



Vision Test Record

F/CCS/10
Revision No. 00
Page 1 of 1

CANDIDATE NAME:..... IDENTIFICATION NO:..... COMPANY:.....

NEAR VISION ACUITY TEST	COLOR CONTRAST DIFFERENTIATION TEST	SHADES OF GRAY TEST (for RT-FI, Radiographic Film Interpretation candidate only)
CRITERIA		
Shall be able to read Jaeger number 1 (J1) at not less than 30.5 cm.	Shall be capable of adequately distinguishing and differentiating colors in in Ishihara color chart.	Shall be able to read gray numbers from Shades of Gray developed by Dr.Kolbl One/TUV/BV
RESULT		
Writing record:.....	Number of Ishihara plates correctly interpreted: <input type="text"/> (MINIMUM OF FIRST 17)	Number of gray numbers correctly read: <input type="text"/> (MINIMUM OF 20 FROM 25)
<input type="checkbox"/> Meets the requirement without correction <input type="checkbox"/> Meets the requirement with correction <input type="checkbox"/> Does not meet the requirement	<input type="checkbox"/> Passes an Ishihara test <input type="checkbox"/> Fails an Ishihara test In case of failure of an Ishihara test, detail of colour which observed as deficiency shall be mentioned here:.....	<input type="checkbox"/> Meets the requirement <input type="checkbox"/> Does not meet the requirement
DETAILS OF PERSON CARRYING OUT AND RECORDING ANY OF THE ABOVE TESTS		
Signature & Stamp: Full Name:..... Appointment/Title:.....	Date of test:/...../..... Expiry date:/...../..... (maximum 12 months from the date of test)	NOTE: (if any)